



Athens 01-08-2013

Determination of lycopene content in raw tomatoes

Plant material: *Solanum lycopersicum*, type "elpida"

Producer of plant material: AGAN EPE

Method of analysis: A spectrophotometric method was employed for the determination of lycopene content. Samples were first chopped and homogenized in a laboratory homogenizer. Approximately 0.3 to 0.6 g of the homogenized samples were weighed in Erlenmayer flasks and 5 mL of 0.05% (w/v) BHT in acetone, 5 mL of ethanol and 10 mL of hexane were added. The flasks were placed in an ice bath and stirred on a magnetic stirring plate for 15 min. After shaking, 3 mL of deionized water were added to each flask and the samples were shaken for 5 min on ice. Samples were then left at room temperature for 5 min to allow the separation of both phases. The absorbance of the hexane layer (upper layer) was measured in a 1-cm-path-length quartz cuvette at 503 nm using hexane as a blank. The lycopene content was calculated using its molar extinction coefficient in hexane ($17.2 \cdot 10^4 \text{ M}^{-1} \cdot \text{cm}^{-1}$) determined in the literature (Markovic, K.; Hruskar, M.; Vahcic, N. *Nutr. Res.* **2006**, 26, 556–560). In this case, the Lambert–Beer law can be described as:

Absorbance at 503 nm (A_{503}) = $\epsilon(\text{M}^{-1} \cdot \text{cm}^{-1}) \cdot b(\text{cm}) \cdot [\text{Lycopene concentration (M)}]$

By properly substituting the molar extinction coefficient of lycopene in hexane ($17.2 \cdot 10^4 \text{ M}^{-1} \cdot \text{cm}^{-1}$), as well as its molecular weight (536.9 g) and by changing the units, the final equation is:

Lycopene content (mg/kg) = $A_{503} \cdot 31.2/\text{g tissue}$

The analysis was performed in five individual tomato samples.

Lycopene content (mg/kg): 49.4 ± 4.7